



## Laboratory Collection Guide for Non-blood Specimens

Refer to our electronic test catalog for our full test menu and complete collection instructions.

<https://unitypointfinley.testcatalog.org>

Test	Collection Device Name	Special Notes
Bacterial Cultures (aerobe and anaerobe)  Fungal Cultures  Strep A PCR	ESwab  	The caps may be white, blue, or green depending on the swab size.  Write the source on the patient label. Example: Left leg wound.
Chlamydia  Gonorrhea  Trichomonas  Vaginitis Molecular testing (≥14 yrs)	GeneXpert Collection Swab  	Use for female genital specimens, male/female rectal specimens, and male/female throat specimens.  For male patients: chlamydia, gonorrhea, and trichomonas is not available on genital specimens. Submit a urine specimen for male patients.

20250613

20250613

À|à\|À| ^ÀFHJÍ G ÈSæ ÒÁçä, ^àÀ} ÀGHEDG ÈÜ|ä çàÀ} ÀGHEDG ÁFKFÀE ÀÒÙVÈÜæ ^Á Á Á